PreScan/ReScan에 방해가 되는 자가형광 대처법 (Applicable Autofluorescence Quenching Methods for PreScan/ReScan)

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Classification	Causes
Endogenous autofluorescence	Flavins (Latin flavus meaning "yellow"): riboflavin, FMN, FAD, FADH/FADH ₂ Porphyrins (RBC hemoglobin originated, 600-730 nm) NAD(P)H (420-550 nm) Lipofuscin (420-700 nm) Elastic fibers (420-600 nm)
Extrinsic fluorescence	Chemical agents, drugs, natural product for in vivo pharmacological uses (emission spectrum varies dependent upon the origin species)
Aldehyde-induced fluorescence	Formaldehyde crosslinking with amino-group seems create broad green to yellow (450-580 nm) fluorescence overall entire tissue.

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Method	Assessment
Time-gated detection method (FLIM and ADOTA fluorophore)	Excellence in separating fluorescent signals Access to the equipment and operation is quite complicated and expensive
Chemical quenching kit (Vector®TrueView™ SP-8400)	Easy to use and more than moderate autofluorescence quenching capability Price per kit (up to 100 tissue slides) is around 200 USD
Sudan Black B (SBB)	Good at 0.1% of SBB in frozen section Moderate at 0.1% of SBB in paraffin section SBB also dissolves many types of lipids on the tissue
LED or laser-induced photobleaching	Prolonged sample preparation time (6-12 hours for photobleaching under 25C°) Delayed exposure time for thicker-sectioned materials

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