


Seahorse XFe96 easy-to-follow guide

**D-2**

Cell seeding  
(4,000 – 8,000 cells/well)



**X** A1, A12, H1, H12  
calibration well only  
(No cell seeding)  
Instead, top off with media

**D-1**

Sensor cartridge hydration  
(Calibrant soaking:  
200µL/well, **37°C/non-CO<sub>2</sub>**  
overnight)

Seahorse XFe96  
**Turn ON** (overnight)

Wave software ON  
Input experimental design

Step  
↔  
exchangeable

**D-day**

Prepare **XF assay media**  
(50mL/plate, **37°C/non-CO<sub>2</sub>**)

Prepare **compound stock**  
solution (A, B, C or A, B)

**Dispense compound** stock  
solution into drug ports (A, B,  
C, D) -> Keep it under  
**37°C/non-CO<sub>2</sub>** (30 min)

Wave software ON  
Input experimental design

**Start Run: sensor cartridge**  
**calibration** (30 min)

**Cell culture plate -> XF assay**  
**media change**

1. **Rinse 2 times** with  
100µL/well of XF assay  
media or 1X PBS
2. **Top off each well** with  
180 µL of XF assay media  
**37°C/non-CO<sub>2</sub>** (40 min)


**Phase contrast imaging with**  
**Cytation1** (20 min)

**Plate switching** (calibration  
plate -> cell culture plate)

Click **"I'm ready"** (1 hr 40 min)

**Normalization** (20 min)

**Data export** (Excel)



Mitochondrial Respiration

