Single cell tracking을 통한 migration ability analysis (Operetta CLS, PerkinElmer)

세포체학실험실에 보유 중인 Operetta CLS (PerkinElmer, MA, US)는 high content analysis (HCA)가 가능한 96-well 또는 384-well 기반 imaging screening 장비로서, well-based & cell-based fluorescence quantification, target cell counting & distribution analysis, apoptosis rough screening & cell scoring without specific fluorescent marker (with Hoechst only), cell morphology & classification, texture (specific superficial features on cells) analysis 및 machine learning을 통한 분석 algorithm building block까지 대량(multi-well plate기반)의 sample의 screening에 활용되고 있습니다. 뿐만 아니라, 실험군과 대조군간에 Z-factor를 자동계산하여 제공하므로써, 해당실험이 추후 대량 scale-up이 가능한지 등의 유효성 여부를 판단하는데 도움을 주고 있습니다. 이번에 살펴 볼 single cell tracking은 HCA장비가 제공하는 분석 tool중에서 cell fluorescent marker 유무에 관계없이 세포의 이동능력을 평가하는 손쉬운 방법입니다.

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
Image Analysis Results												
Results		Summary		Properties Cells								
Field	Number of Tracked Cells	Age[s]	Current Step [μm]	Current Speed [$\mu\text{m/s}$]	Mean Square Displacement	Current Displacement X [μm]	Current Displacement Y [μm]	Number of Tracks	Number of Timepoints per Track	Accumulated Distance [μm] per Track	Displacement [μm] per Track	Average Speed [$\mu\text{m/s}$] per Track
1	179	4863.02	6.17968	0.0171666	4492.59	4.29472	8.06232	475	9.55789	52.0547	35.1853	0.0175587

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