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

분석 algorithm:

Open file (IF형광 image) -> Separate RGB -> Adjustment of inclusive threshold of "DAPI image" by input of low threshold -> Move to the next tap "Clip" and check "Create binary mask", "1-bit depth" and "Apply" -> 1bit black & white image가 생성 됨.


Arithmetic:

Go to "Arithmetic" -> Source image 1: 세포(cell) image 선택 -> Source image 2: 1bit black & white image 선택 -> Image mode: AND -> Bit depth: 16bit -> Operation: "Logical AND (핵 부위 오려내기) or Subtract (세포질 부위 오려내기)" -> Apply

Threshold Adjustment:

Go to "Threshold Image" -> Select source image by clicking "Cytosol" or "Nuclei" image -> Inclusive threshold -> Decide "Low" threshold by controlling a sliding bar. Tip: in order to make reliable comparison of batch images, input a fixed number onto the Low threshold box in inclusive threshold mode.

Cell_Green intensity			
	Area	Average intensity	Total intensity
Average	441.475	58.87	35795.125
Std. Dev.	12835	6.245	1032934.313
Minimum	2	50	100
Maximum	410724	133.538	33048396
Total	452070	60283.148	36654208

Nucleus_Green intensity				N/C ratio
	Area	Average intensity	Total intensity	%
Average	632.692	95.93	62475.703	
Std. Dev.	1308.21	39.028	131943.281	
Minimum	2	50.5	101	
Maximum	13895	241	1373481	
Total	205625	31177.412	20304604	55.395

Integrated Morphometry Analysis:

Go to IMA (Integrated Morphometry Analysis) -> Select source image -> Check segment using mask and select the image for masking (Masking is a useful tool to extract specific region for its analysis. If you mask your image with DAPI, only the image region which covered by DAPI signal will be analyzed) -> Decide parameters desired (check display boxes) -> For more specific selections, check "Filter" and give Limit1 and Limit2 numbers. -> Measure -> Confirm object data and summary -> Data export.