

FACSAria II Detector Configuration

2020.11.20

Laser	Dectector Name	BP Filter	LP Filter	Fluorochrome Detected
Blue laser (488nm)	A	780/60	755LP	PE-Cy7
	B	710/50	685LP	PerCP-Cy5.5, BB700
	C	660/20	635LP	PE-Cy5
	D	610/20	600LP	PE-Texas Red
	E	575/25	550LP	PE
	F	530/30	505LP	FITC, BB515, Alexa 488, GFP
	G	488/10		SSC
Red laser (640nm)	A	780/60	750LP	APC-Cy7, APC-H7
	B	730/45	690LP	Alexa Fluor 700, APC-R700
	C	670/30		APC
Violet laser (405nm)	A	780/60	750LP	BV786
	B	710/50	690LP	BV711
	C	670/30	630LP	BV650
	D	610/20	595LP	BV605
	E	525/50	505LP	BV510
	F	450/40		BV421